

Product Datasheet

Normal Rabbit IgG Control antibody (orb633278)

Description	Rabbit polyclonal antibody to Normal Rabbit IgG Control
Conjugation	Unconjugated
Tested Applications	Control
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, pH7.4. 4°C for 24 months.
Note	For research use only
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Expiration Date	12 months from date of receipt.