

Product Datasheet

Affinipure Rabbit Anti-human Fibrinogen antibody (Cy3) (orb633280)

Description	Rabbit polyclonal antibody against Affinipure Rabbit Anti-human Fibrinogen conjugated to Cy3
Conjugation	Cy3
Tested Applications	FC, IF
Preservatives	5mg/ml rAlbumin (IgG-Free, Protease-Free)
Form/Appearance	liquid
Storage	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.
Note	For research use only
Hazard Information	This product is for research use only, not for diagnostic or therapeutic use.
Expiration Date	12 months from date of receipt.