



Product Datasheet

LTB4-R1 antibody (orb6335)


Description	LTB4-R1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human LTB4-R1 (101-200/352 aa)
Target	LTB4-R1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q15722
Expiration Date	12 months from date of receipt.



Paraffin-embedded mouse lung tissue fixe...



Paraffin-embedded human endometrium carc...



Western Blotting analysis of mouse heart...