

Product Datasheet

CHRNA6 Antibody (orb636116)

Description	Mouse monoclonal antibody to CHRNA6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC
Immunogen	Purified recombinant fragment of human CHRNA6 (AA: 26-239) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	57kDa
Entrez	8973
Dilution Range	IHC: 1/200 - 1/1000; ICC: 1/100 - 1/500; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.