

Product Datasheet

AP2M1 Antibody (orb636133)

Description	Mouse monoclonal antibody to AP2M1
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified recombinant fragment of human AP2M1 (AA: 298-435) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	49.7kDa
Entrez	1173
Dilution Range	WB: 1/300 - 1/800; ELISA: 1/10000
Expiration Date	12 months from date of receipt.