

Product Datasheet

NRP1 Antibody (orb636134)

Description	Mouse monoclonal antibody to NRP1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified recombinant fragment of human NRP1 (AA: 22-644) expressed in HEK293 cells.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	103.1kDa
Entrez	8829
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.