

Product Datasheet

ADORA2A Antibody (orb636188)

Description	Mouse monoclonal antibody to ADORA2A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC
Immunogen	Purified recombinant fragment of human ADORA2A (AA: 274-412) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	44.7kDa
Entrez	135
Dilution Range	ICC: 1/100 - 1/500; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.