



Product Datasheet

GRM2 Antibody (orb636229)


Description	Mouse monoclonal antibody to GRM2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	Purified recombinant fragment of human GRM2 (AA: extra 414-558) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Form/Appearance	Purified antibody in PBS with 0.05% sodium azide.
Concentration	1mg/ml
Storage	This product will ship in a box containing blue ice at a temperature of 4°C.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	95.6kDa
Entrez	2912
Dilution Range	ELISA: 1/10000, FACS: 1/200 - 1/400, WB: 1/500 - 1/2000
Expiration Date	12 months from date of receipt.



Western blot analysis of SK-N-SH (1), Ju...



Line graph illustrates about the Ag-Ab r...



Flow cytometric analysis of SK-N-SH cell...