

Product Datasheet

EDA-m7GTP (ATTO-540Q) (orb63636)

Description	EDA-m7GTP conjugated to ATTO-540Q.
Form/Appearance	solution in water; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 543 nm, ϵ 105.0 L mmol ⁻¹ cm ⁻¹ , (Tris-HCl pH 7.5).
Formula	C ₁₄ H ₂₄ N ₇ O ₁₅ P ₃ -ATTO 540Q
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1163.30 g/mol (free acid)
SMILES	[O-]]P(=O)(O)OP(=O)(O)OP(=O)(O)OC[C@H]1O[C@@H](n2c[n+](c3c(=O)[nH]c(nc23)N)C)[C@H](O)[C@@H]1OC(=O)NCCN*.[O-]]P(=O)(O)OP(=O)(O)OP(=O)(O)OC[C@H]1O[C@@H](n2c[n+](c3c(=O)[nH]c(nc23)N)C)[C@H](OC(=O)NCCN*)[C@@H]1O
Expiration Date	12 months from date of receipt.