

Product Datasheet

TRIM29 Antibody (orb636395)

Description	Mouse monoclonal antibody to TRIM29
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Purified recombinant fragment of human TRIM29 (AA: 451-588) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	65.8kDa
Entrez	23650
Dilution Range	WB: 1/500 - 1/2000; ICC: 1/100 - 1/500; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.