



Product Datasheet

ULBP2 Antibody (orb636631)


Description	This gene encodes a major histocompatibility complex (MHC) class I-related molecule that binds to th...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Purified recombinant fragment of human ULBP2 (AA: 26-216) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	Expected MW is 47.6 kDa
Entrez	80328
Dilution Range	WB: 1/500 - 1/2000 ICC (Immunocytochemistry): 1/500 - 1/2000 FC: 1/200 - 1/400 ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Western blot analysis using ULBP2 mAb ag...



Western blot analysis using ULBP2 mAb ag...



Black line:
Control
Antigen (100
ng);Pur...