



Product Datasheet

CD364 Antibody (orb636638)


Description	Mouse monoclonal antibody to CD364
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	Purified recombinant fragment of human CD364 (AA: 79-210) expressed in E. Coli.
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	49.5kDa
Entrez	221476
Dilution Range	WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
HL-60 cells ...



Western blot
analysis of
HEK293 (1)
and ...



Western blot
analysis of
THP-1 (1),
PC-3...