

Product Datasheet

Human PSCA Protein (orb636720)

Description

Human PSCA, His Tag (PSA-H52H6) is

contains AA Leu 12 - ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

PSCA

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 10.4 kDa. The protein migrates as 18-25 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_005663.2

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human PSCA, His Tag (orb636720) is expressed from human 293 cells (HEK293). It contains AA Leu 12 - Ser 86 (Accession # O43653-1).

MW

10.4 kDa

Uniprot ID

O43653

NCBI

NP_005663.2

Expiration Date

6 months from date of receipt.