



Product Datasheet

ELNE antibody (orb636861)


Description	Rabbit polyclonal antibody to ELNE
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IF: 1: 100-1:500, ICC: 1: 100-1:500, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cells...



IHC-P image of rat spleen tissue ELNE antibody...



Immunohistochemical staining of paraffin-embedded tissue...