

Product Datasheet

RAB27B antibody (orb637040)

Description	Rabbit polyclonal antibody to RAB27B
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:100-1:500, ICC: 1:100-1:500
Expiration Date	12 months from date of receipt.