

Product Datasheet

EDA-GTP gamma S (ATTO-700) (orb63715)

Description	EDA-GTP gamma S conjugated to ATTO-700.
Form/Appearance	solution in water; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 700 nm, λ_{em} 716 nm, ϵ 120.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₄₃ H ₅₅ N ₁₀ O ₁₉ P ₃ S ₂
Purity	≥ 90% (HPLC)
MW	Theoretical MW: 1173.01 g/mol (free acid); Detected MW: 1172.23 g/mol (free acid)
SMILES	<chem>N(CCNC(=O)O[C@@H]1[C@@H](COP(=O)(O)OP(=O)(O)OP(=S)(O)O)[C@@H](n2cnc3c(=O)[nH]c(nc23)N)[C@@H]1O)C(=O)CCCN1c2cc3Oc4cc5=[N+](CC)C(C=C(c5cc4=Nc3cc2C(=CC1(C)C)CS(=O)(=O)[O-])C)(C)C.N(CCNC(=O)O[C@@H]1[C@@H](n2cnc3c(=O)[nH]c(nc23)N)O[C@@H](COP(=O)(O)OP(=O)(O)OP(=S)(O)O)[C@@H]1O)C(=O)CCCN1c2cc3Oc4cc5=[N+](CC)C(C=C(c5cc4=Nc3cc2C(=CC1(C)C)CS(=O)(=O)[O-])C)(C)C</chem>
Expiration Date	12 months from date of receipt.