



Product Datasheet

HRP2 antibody (orb637262)


Description	Rabbit polyclonal antibody to HRP2
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from H...



Western blot analysis of Lane 1: VERO tr...



Immunohistochemical staining of Rat test...