

## Product Datasheet

### CRISP3 antibody (orb637894)

|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CRISP3                  |
| <b>Species/Host</b>        | Rabbit                                                |
| <b>Reactivity</b>          | Human, Monkey, Mouse                                  |
| <b>Conjugation</b>         | Unconjugated                                          |
| <b>Tested Applications</b> | ICC, IF, WB                                           |
| <b>Note</b>                | For research use only                                 |
| <b>Clonality</b>           | Polyclonal                                            |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IF: 1: 100-1:500, ICC: 1: 100-1:500 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                       |