

Product Datasheet

MAP4K4 (phospho-Ser631) antibody (orb6380)

Description	MAP4K4 (phospho-Ser631) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MAP4K4 around the phosphorylation site of Ser631 (SR(p-S)PV)
Target	phospho-MAP4K4 (Ser631)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	142 kDa
Uniprot ID	O95819
Expiration Date	12 months from date of receipt.