

Product Datasheet

EDA-ATP gamma S (ATTO-665) (orb63808)

Description	EDA-ATP gamma S conjugated to ATTO-665.
Form/Appearance	solution in water; Color: blue; pH: 7.5 ± 0.5
Concentration	0.10 mM-0.11 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ _{exc} 662 nm, λ _{em} 680 nm, ε 160.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₁₃ H ₂₂ N ₇ O ₁₃ P ₃ S-ATTO 665
Purity	≥ 90% (HPLC)
MW	Theoretical MW: 1213.33 g/mol (free acid)
SMILES	N(CCNC(=O)O[C@@H]1[C@@H](COP(=O)(O)OP(=O)(O)OP(=S)(O)O)[C@@H](n2cnc3c(N)ncnc23)[C@@H]1O)*.N(CCNC(=O)O[C@H]1[C@H](n2cnc3c(N)ncnc23)O[C@H](COP(=O)(O)OP(=O)(O)OP(=S)(O)O)[C@H]1O)*
Expiration Date	12 months from date of receipt.