

Product Datasheet

Nicotinic Acetylcholine Receptor alpha 7 antibody (orb638225)

Description	Rabbit polyclonal antibody to Nicotinic Acetylcholine Receptor alpha 7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:500-1:2000, IF: 1: 100-1:500, ICC: 1: 100-1:500, IHC-P: 1:50-1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.