

Product Datasheet

GNDF Receptor alpha 2 antibody (orb638285)

Description	Rabbit polyclonal antibody to GDNF Receptor alpha 2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	IF: 1: 100-1:500, ICC: 1: 100-1:500, WB: 1:500-1:2000, IHC-P: 1:50-1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.