

## Product Datasheet

### Synaptojanin antibody (orb638544)

|                            |                                                                               |
|----------------------------|-------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Synaptojanin                                    |
| <b>Species/Host</b>        | Rabbit                                                                        |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                  |
| <b>Tested Applications</b> | ELISA, ICC, IF, WB                                                            |
| <b>Note</b>                | For research use only                                                         |
| <b>Clonality</b>           | Polyclonal                                                                    |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IF: 1: 100-1:500, ICC: 1: 100-1:500, ELISA: 1:20000-1:40000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                               |