

## Product Datasheet

### PITX2 antibody (orb638552)

|                            |                                          |
|----------------------------|------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PITX2      |
| <b>Species/Host</b>        | Rabbit                                   |
| <b>Reactivity</b>          | Human, Mouse, Rat                        |
| <b>Conjugation</b>         | Unconjugated                             |
| <b>Tested Applications</b> | ELISA, WB                                |
| <b>Note</b>                | For research use only                    |
| <b>Clonality</b>           | Polyclonal                               |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, ELISA: 1:20000-1:40000 |
| <b>Expiration Date</b>     | 12 months from date of receipt.          |