

Product Datasheet

C-MYC-tag antibody (orb638732)

Description	Rabbit polyclonal antibody to C-MYC-tag
Species/Host	Rabbit
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB: 1:3000-1: 10000, IF: 1:200, ICC: 1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.