

Product Datasheet

mCherry-Tag antibody (orb638738)

Description	Mouse monoclonal antibody to mCherry-Tag
Species/Host	Mouse
Reactivity	All, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, WB
Note	For research use only
Clonality	Monoclonal
Dilution Range	WB: 1:5000, IF: 1:200, ICC: 1:200, ELISA: 1:20000-1:40000
Expiration Date	12 months from date of receipt.