



Product Datasheet

SARS N (nucleocapsid) antibody (orb638864)

Description	SARS N (nucleocapsid) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant SARS N (nucleocapsid) (1-422/422aa)
Target	SARS-CoV Nucleocapsid
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	≥1mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	P59594
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Lane 1: SARS-CoV-2 N



Western blot analysis of Lane 1: 2019-nCoV N