


Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid protein (orb638879)

Description	Nucleocapsid (N) protein is the most abundant protein found in coronavirus. CoV N protein is a highl...
Reactivity	Virus
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Nucleocapsid protein
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 47.3 kDa. The protein migrates as 60-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	QHO62115.1
Purity	95%
Source	SARS-CoV-2 Nucleocapsid protein, His Tag (orb638879) is expressed and secreted by human 293 cells (HEK293). It contains AA Met 1 - Ala 419 (Accession # QHO62115.1).
MW	47.3 kDa
NCBI	QHO62115.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of 2019-nCoV (COVID-19...