

## Product Datasheet

### COVID-19 S1 CTD protein (orb638895)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Recombinant COVID-19 S1 Protein CTD is produced by mammalian expression system with a 8His; Flag tag...  |
| <b>Reactivity</b>          | Virus  |
| <b>Endotoxins</b>          | Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.   |
| <b>Tested Applications</b> | ELISA, MS, SDS-PAGE, WB  |
| <b>Preservatives</b>       | Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.   |
| <b>Form/Appearance</b>     | Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.   |
| <b>Storage</b>             | Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles. |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Please minimize freeze-thaw cycles.  |
| <b>Purity</b>              | Greater than 95% as determined by reducing SDS-PAGE.   |
| <b>Source</b>              | Human cells  |
| <b>Entrez</b>              | <a href="#">43416, 1</a>   |
| <b>Expiration Date</b>     | 6 months from date of receipt.   |