

Product Datasheet

TEAD1 antibody (orb638927)

Description	Rabbit polyclonal antibody to TEAD1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human TEAD1(NP_068780.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50KDa
Uniprot ID	P28347
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF: 1:50-200, ChIP: 1:20-50
Expiration Date	12 months from date of receipt.