

Product Datasheet

Beclin 1 3D antibody (orb638951)

Description	Rabbit polyclonal antibody to Beclin 1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Beclin 1(NP_003757.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60KDa
Uniprot ID	Q14457
Dilution Range	WB: 1:500-2000, IF: 1:50-200
Expiration Date	12 months from date of receipt.