

Product Datasheet

EDA-AppNHp (ATTO-740) (orb63898)

Description	EDA-AppNHp conjugated to ATTO-740.
Form/Appearance	solution in water; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 743 nm, λ_{em} 763 nm, ϵ 120.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₁₃ H ₂₃ N ₈ O ₁₃ P ₃ -ATTO 740
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1041.29 g/mol (free acid)
SMILES	[C@@H]1([C@H]([C@@H](O[C@@H]1COP(=O)(OP(=O)(NP(=O)(O)O)O)O)n1c2c(nc1)c(ncn2)N)OC(=O)NCCN*)O.[C@@H]1([C@H]([C@@H](O[C@@H]1COP(=O)(OP(=O)(NP(=O)(O)O)O)O)n1c2c(nc1)c(ncn2)N)O)OC(=O)NCCN*
Expiration Date	12 months from date of receipt.