

Product Datasheet

ARHGAP4 antibody (orb639025)

Description	Rabbit polyclonal antibody to ARHGAP4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human ARHGAP4(NP_001158213.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105kDa
Uniprot ID	P98171
Dilution Range	WB: 1:500-2000, IHC-P: 1:50-200, IF: 1:50-200
Expiration Date	12 months from date of receipt.