

Product Datasheet

EDA-AppNHp (ATTO-633) (orb63905)

Description	EDA-AppNHp conjugated to ATTO-633.
Form/Appearance	solution in water; Color: blue; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ _{exc} 630 nm, λ _{em} 651 nm, ε 130.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₄₈ H ₆₂ N ₁₁ O ₁₅ P ₃
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1126.01 g/mol (free acid); Detected MW: 1125.36 g/mol (free acid)
SMILES	<chem>N(CCNC(=O)O[C@@H]1[C@@H](COP(=O)(O)OP(=O)(O)NP(=O)([O-])O)O[C@@H](n2cnc3c(N)ncnc23)[C@@H]1O)C(=O)CCCN(C(=O)c1cccc1C1=C2C=CC(=[N+](C)C)C=C2C(c2c1cc1c(c2)N(CC)CCC1)(C)C)C.O=C(O[C@@H]1[C@@H](n2cnc3c(N)ncnc23)O[C@@H](COP(=O)(O)OP(=O)(O)NP(=O)([O-])O)[C@@H]1O)NCCNC(=O)CCCN(C(=O)c1cccc1C1=C2C=CC(=[N+](C)C)C=C2C(c2c1cc1c(c2)N(CC)CCC1)(C)C)C</chem>
Expiration Date	12 months from date of receipt.