

Product Datasheet

ITSN1 antibody (orb639060)

Description	Rabbit polyclonal antibody to ITSN1
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human ITSN1(NP_001317941.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	137kDa
Uniprot ID	Q15811
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.