

Product Datasheet

SIPA1L1 antibody (orb639126)

Description	Rabbit polyclonal antibody to SIPA1L1
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SIPA1L1(NP_056371.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	180kDa
Uniprot ID	O43166
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.