

Product Datasheet

BTBD10 antibody (orb639187)

Description	Rabbit polyclonal antibody to BTBD10
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human BTBD10.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9BSF8
Entrez	84280
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.