

Product Datasheet

FBF1 antibody (orb639194)

Description	Rabbit polyclonal antibody to FBF1
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human FBF1(Q8TES7).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54kDa/69kDa/98kDa/125kDa/127kDa
Uniprot ID	Q8TES7
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.