

Product Datasheet

USH1G antibody (orb639213)

Description	Rabbit polyclonal antibody to USH1G
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human USH1G(NP_775748.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q495M9
Dilution Range	WB: 1:500-2000, IF: 1:50-200
Expiration Date	12 months from date of receipt.