

## Product Datasheet

### PUSL1 antibody (orb639214)

|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PUSL1                     |
| <b>Species/Host</b>        | Rabbit                                                  |
| <b>Reactivity</b>          | Human, Rat                                              |
| <b>Conjugation</b>         | Unconjugated                                            |
| <b>Tested Applications</b> | WB                                                      |
| <b>Immunogen</b>           | Recombinant protein of human PUSL1.                     |
| <b>Preservatives</b>       | Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.             |
| <b>Note</b>                | For research use only                                   |
| <b>Isotype</b>             | IgG                                                     |
| <b>Clonality</b>           | Polyclonal                                              |
| <b>Purity</b>              | Affinity purification                                   |
| <b>Uniprot ID</b>          | <b>Q8N0Z8</b>                                           |
| <b>Entrez</b>              | <b>126789</b>                                           |
| <b>Dilution Range</b>      | WB: 1:500-2000                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                         |