

## Product Datasheet

### PLEKHA7 antibody (orb639223)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PLEKHA7                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                         |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Tested Applications</b> | ICC, IF, WB                                               |
| <b>Immunogen</b>           | Recombinant fusion protein of human PLEKHA7(NP_778228.3). |
| <b>Preservatives</b>       | PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.          |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.               |
| <b>Note</b>                | For research use only                                     |
| <b>Isotype</b>             | IgG                                                       |
| <b>Clonality</b>           | Polyclonal                                                |
| <b>Uniprot ID</b>          | <b>Q6IQ23</b>                                             |
| <b>Dilution Range</b>      | WB: 1:500-2000, IF: 1:50-200                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |