

Product Datasheet

C14orf37 antibody (orb639224)

Description	Rabbit polyclonal antibody to C14orf37
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human C14orf37(NP_001001872.2).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86TY3
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.