

Product Datasheet

TUSC1 antibody (orb639241)

Description	Rabbit polyclonal antibody to TUSC1
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human TUSC1(NP_001004125.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q2TAM9
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.