

Product Datasheet

TAZ1 antibody (orb639245)

Description	Rabbit polyclonal antibody to TAZ1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of mouse TAZ1(NP_001186665.1).
Preservatives	PBS with 0.01% thiomersal,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	191KDa
Uniprot ID	Q8N9V7
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.