

Product Datasheet

POU3F2 antibody (orb639252)

Description	Rabbit polyclonal antibody to POU3F2
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human POU3F2(NP_005595.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55KDa
Uniprot ID	P20265
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.