

Product Datasheet

APOBEC3H antibody (orb639301)

Description	Rabbit polyclonal antibody to APOBEC3H
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human APOBEC3H.
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q6NTF7
Entrez	164668
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.