

Product Datasheet

RTN4R antibody (orb639351)

Description	Rabbit polyclonal antibody to RTN4R
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human NOGO Receptor 1 (RTN4R)(NP_075380.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	66KDa
Uniprot ID	Q9BZR6
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.