

Product Datasheet

SRRT antibody (orb639372)

Description	Rabbit polyclonal antibody to SRRT
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant fusion protein of human SRRT(NP_001122325.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	96kDa/100kDa
Uniprot ID	Q9BXP5
Dilution Range	IHC-P: 1:50-200
Expiration Date	12 months from date of receipt.