

Product Datasheet

GRIN1 (Phospho-S896) antibody (orb639376)

Description	Rabbit polyclonal antibody to GRIN1 (Phospho-S896)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S896 of human NMDAR1
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	120kDa
Uniprot ID	Q05586
Entrez	2902
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.