

Product Datasheet

mTOR (Phospho-S2481) antibody (orb639378)

Description	Rabbit polyclonal antibody to mTOR (Phospho-S2481)
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S2481 of human mTOR(NP_004949.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
MW	289kDa
Uniprot ID	P42345
Dilution Range	WB: 1:500-2000
Expiration Date	12 months from date of receipt.